


INERTIA® ONE-SHOT™ PEDICLE SCREW

A New Direction With Momentum



Probe tip aids
in tracking
the pedicle.



One-Step Delivery

Probe tip allows the surgeon the choice of placing the screw in One-Shot or using an awl, drill and tap to place the screw.

"I wanted a screw that would allow me M.I.S. placement without the need of a cannulated system, targeting needles, k-wires, etc., and after using the Inertia screw I knew I had found a platform to build on.

The One-Shot allows me to place and direct the screw as needed minimizing both my time of targeting the pedicle and risk of error."

— Dennis L. Abernathie, MD



The One-Shot is a great option to use as a revision screw with 5.5, 6.5, 6.75, 7.0 and 7.5mm screw sizes.

NEXXT
Spine

Nexxt Spine
14425 Bergen Blvd.
Suite B
Noblesville, IN 46060

www.NexxtSpine.com
info@NexxtSpine.com
phone: 317.436.7801
fax: 317.245.2518

70-006B